



US008784021B2

(12) **United States Patent**
Luksch et al.

(10) **Patent No.:** **US 8,784,021 B2**
(45) **Date of Patent:** **Jul. 22, 2014**

(54) **JIG DEVICE AND APPARATUS AND METHOD OF MAKING A DENTAL PROSTHESIS OR PATTERN THEREFOR**

(75) Inventors: **Derrick Luksch**, South Beach, OR (US); **Christophorus Luksch**, Corvallis, OR (US)

(73) Assignee: **O'Brien Dental Lab, Inc.**, Corvallis, OR (US)

(*) Notice: Subject to any disclaimer, the term of this patent is extended or adjusted under 35 U.S.C. 154(b) by 1294 days.

(21) Appl. No.: **12/612,872**

(22) Filed: **Nov. 5, 2009**

(65) **Prior Publication Data**

US 2011/0104642 A1 May 5, 2011

(51) **Int. Cl.**
A61C 5/10 (2006.01)
B23C 3/00 (2006.01)

(52) **U.S. Cl.**
USPC **409/132**; 409/164; 409/174; 409/198;
409/221; 433/51; 433/223; 269/7

(58) **Field of Classification Search**
USPC 409/131-132, 164-165, 174, 189,
409/197-198, 219, 221; 433/51, 223; 269/7
IPC A61C 5/10
See application file for complete search history.

(56) **References Cited**

U.S. PATENT DOCUMENTS

2,633,776 A * 4/1953 Schenk 409/132
2,676,499 A * 4/1954 Walther 408/1 R
4,730,382 A * 3/1988 Parsons 29/563
4,904,348 A * 2/1990 Domes et al. 205/67
5,079,821 A * 1/1992 Parsons 29/424

5,116,174 A * 5/1992 Fried et al. 409/79
5,634,250 A * 6/1997 Mihailovic 29/26 A
5,836,729 A * 11/1998 Fink 409/132
6,019,554 A * 2/2000 Hong 409/132
6,082,938 A * 7/2000 Fink 409/132
6,739,959 B2 * 5/2004 Bodenmiller et al. 451/364
6,915,178 B2 7/2005 O'Brien et al.
7,094,012 B1 * 8/2006 Bang et al. 409/221
7,234,938 B2 * 6/2007 Bodenmiller 433/51
7,252,466 B1 * 8/2007 Roch 409/187

FOREIGN PATENT DOCUMENTS

SU 1227373 A * 4/1986

* cited by examiner

Primary Examiner — Daniel Howell

(74) *Attorney, Agent, or Firm* — John J. Connors; Connors & Assoc. PC

(57) **ABSTRACT**

A jig device holding a work piece is used with a milling machine that cuts from the work piece a prosthesis having a surface essentially free of defects. The jig device is inserted repeatedly into the machine in the same predetermined location, even upon inverting the jig device. This maintains the spatial relationship between the jig device and milling so, at start up, they are always in the same exact spatial relationship. The milling machine partially cuts through a first side of the work piece to form in the work piece a cavity having a surface part corresponding to at least a portion of the prosthesis being made. The jig device is then removed from the machine, inverted and reinserted into the machine at the same predetermined location. A support material is placed in the cavity prior to reinserting the jig device. Upon reinserting the jig device, the machine precisely cuts through an opposed side of the work piece to form within the work piece a prosthesis having a surface essentially free of defects. The prosthesis is supported within the work piece by the support material. The prosthesis and support material are separated to complete our method.

10 Claims, 10 Drawing Sheets

